SYSTÈME DE STENT CORONAIRE TITAN OPTIMAX Notice d'utilisation - Français	2
TITAN OPTIMAX CORONARY STENT SYSTEM User manual - English	4
TITAN OPTIMAX KORONARSTENTSYSTEM Benutzungsanleitung - Deutsch	6
SISTEMA STENT CORONARICO TITAN OPTIMAX Istruzioni per l'uso - Italiano	9
SISTEMA DE IMPLANTE CORONARIO TITAN OPTIMAX Instrucciones de uso - Español	11
SISTEMA DE STENT CORONÁRIO TITAN OPTIMAX Manual de utilização - Português	13
TITAN OPTIMAX 冠状动脉支架系统 使用說明書 - 中國	15
ΣΥΣΤΗΜΑ ΣΤΕΦΑΝΙΑΙΟΥ ΣΤΕΝΤ ΤΙΤΑΝ OPTIMAX Οδηγίες χρήσης - Ελληνικά	17
SYSTEM STENTÓW WEWNĘTRZNYCH TITAN OPTIMAX Instrukcja obsługi - Polski	20
KORONÁRNÍ STENT TITAN OPTIMAX Návod k použití - Český	22
SYSTÉM KORONÁRNEHO STENTU TITAN OPTIMAX Informácie pre používateľ'ov - Slovenský	24
TITAN OPTIMAX KORONART STENTSYSTEM Brugervejledning - Dansk	26
TITAN OPTIMAX CORONARY STENT SYSTEM Bruksanvisning - Norsk	28
TITAN OPTIMAX KORONARSTENTSYSTEM Bruksanvisning - Svenska	30
STENT MẠCH VÀNH TITAN OPTIMAX HƯỚNG DẪN SỬ DỤNG - Tiếng Việt	32
TITAN OPTIMAX KOSZORÚÉR- SZTENTRENDSZER Használati útmutató - Magyar	34

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- (1) Tonino P et al., Titanium-Nitride-Oxide-Coated versus Everolimus-Eluting Stents in acute coronary syndrome the randomized TIDES-ACS-Trial, JACC Cardiovascular Intervention
- (2) López-Minguez JR, et al, TIOMAX, A Spanish Multicenter Registry of the real-world use of the Titanium Optimax biostent, Catheter Cardiovascular Intervention
- (3) Karjalainen PP, et al, Clinical outcome of titanium-nitride-oxide-coated cobalt chromium stents in patients with de novo coronary lesions: 12-month results of the OPTIMAX first-in-man study, Catheter Cardiovascular Intervention

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**READ ALL INSTRUCTIONS CAREFULLY BEFORE USE.
FOLLOW ALL INSTRUCTIONS AND PRECAUTIONS IN THE INSTRUCTIONS
TO AVOID ANY POSSIBLE COMPLICATIONS.**

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Sterile. Sterilised with Ethylene Oxide. Pyrogen-free.
Single use. Do not re-sterilise, reprocess and/or re-use.
Do not use if package has been damaged or opened.
Use before the expiry date.
After use, dispose of product and packaging in accordance with local regulations.

TITAN OPTIMAX coronary stent system manual.
©2011 HEXACATH - All rights reserved
TITAN OPTIMAX, Tinox, TiNO, Hexacath and their respective
logos are registered trademarks of Hexacath S.A.S

-I- DESCRIPTION

The TITAN OPTIMAX is a coronary pre-mounted stent system (stent crimped on the balloon of carrier catheter).

The carrier catheter introduces and deploys the stent into the native coronary arteries or vessel grafts of the targeted patients. The TITAN OPTIMAX stent is intended to be implanted in the native coronary arteries or vessel grafts of targeted patients for improving coronary luminal diameters.

The stent structure, made of L605 cobalt-chromium (an alloy of cobalt*, chromium, tungsten and nickel) is a flexible tubulo-modular structure entirely coated with Titanium Nitride Oxide and deployable by means of a balloon catheter. The surface area in contact with the artery varies from 9 mm² to 69 mm² depending on the diameter and length of the stent

The delivery catheter on which the stent is pre-mounted is mainly composed of a flexible distal tip, balloon and tubes made from polymers (Polyamid, Polyether) and a stainless steel hypotube. The outer surface is coated with a non-active hydrophilic coating.

The balloon has two radiopaque markers made from platinum/iridium, proximal and distal, to identify the location of the stent on the balloon. Two proximal markers located at 90 and 100 cm from the distal end help to assess the position of the catheter in relation to the tip of a guide catheter used for radial, brachial or femoral approaches.

The delivery catheter is equipped with a proximal luer-lock port made from polycarbonate for connection to an inflation device.

*Note: the presence of Cobalt material representing more than 0.1% of the stent structure weight classifies it as CMR (carcinogenic, mutagenic or toxic to reproduction) type 1B according to regulation (CE) 1272/2008. However this substance is unlikely to come in contact with the human body, or to be released in the body tissues due to the existence of the Titanium Nitride Oxide coating which covers the entire surface of the stent. Therefore, no particular precaution has to be taken due to the presence of the Cobalt.

-II- OVERVIEW AND REFERENCES

Packaged in a peelable pouch (sterile barrier system) and a cardboard case.

Contents: a TITAN OPTIMAX stent pre-mounted on a catheter, an instruction leaflet, a compliance card, an implant card to be given to the patient.

The TITAN OPTIMAX stent is available in the following lengths: 7, 10, 13, 16, 19, 22, 25, 28, 32, 38 mm and with the following diameters: 2.0, 2.25, 2.5, 2.75, 3.0, 3.5, 4.0, 4.5 and 5.0 mm.

-III- INDICATIONS

The TITAN OPTIMAX Coronary Stent System is indicated for improving coronary luminal diameters in adult patients with symptomatic ischemic heart disease including patients with acute myocardial infarction (STEMI/NSTEMI, unstable angina) and patients with concomitant diabetes mellitus due to de novo coronary artery lesions. The treated lesion length should be less than the nominal stent length (7mm to 38 mm) with a reference vessel diameter of 2.0 to 5.0 mm.

-IV- CONTRAINDICATIONS

- Patients with contraindications to anticoagulant/ antiplatelet therapy
- Allergy to contrast agents.
- Patient with a known hypersensitivity or allergy to Titanium Nitride Oxide (TiNO).
- Major surgery within the previous two weeks, childbirth, puncture of an incompressible vessel or biopsy.
- History of bleeding.
- Pregnancy.
- Lesions proximal to an untreatable segment, significantly limiting blood flow.
- Fibrous or calcified lesions, refractory to dilation (resistant to pre-dilatations at high pressure of 20 Bars).
- In stent restenosis in coronary artery lesions.
- Implantation in tortuous vessels with limited access to the site to be treated, diffuse vascular disease with reduced flow, and lesions located after acute angles, can be challenging situations for stenting.

-V- ADVERSE EFFECTS**5.1 Observed Adverse Effects**

Observed adverse effects comes from the TIDES-ACS clinical trial, the TIOMAX registry and the OPTIMAX "first in man" study. They are summarized in the table below :

	TIDES-ACS	TIOMAX	OPTIMAX first in man
1 month follow-up			
Death		1.4% (7/511)	
Cardiac death		1% (5/511)	
Stent thrombosis (definite or probable)		0.4% (2/511)	
12 months follow-up			
All-cause Death	0.9% (9/989)	4.1% (21/511)	2.2% (5/224)
Non cardiac Death	0.4% (4/989)	2.3% (12/511)	0.9% (2/224)
Cardiac Death	0.5% (5/989)	1.8% (9/511)	1.3% (3/224)
Stent thrombosis (definite or probable)	1.1% (11/989)	0.6% (3/511)	0% (0/224)
Myocardial infarction (MI)	1.8% (18/989)	2.5% (13/511)	3.1% (7/224)
Target lesion revascularisation (TLR)	5.4% (53/989)	2.9% (15/511)	3.1% (7/224)
Target vessel revascularisation (TVR)	6.6% (65/989)	2.9% (15/511)	4.0% (9/224)
Major bleeding	1.2% (12/989)		
18 months follow-up			
All-cause Death	1.1% (11/989)		
Non cardiac Death	0.5% (5/989)		
Cardiac Death	0.6% (6/989)		
Stent thrombosis (definite or probable)	1.1% (11/989)		
Myocardial infarction (MI)	2.2% (22/989)		
Target lesion revascularisation (TLR)	5.8% (57/989)		
Target vessel revascularisation (TVR)	7.1% (70/989)		
Major bleeding	1.4% (14/989)		

- (1) Tonino P, et al., Titanium-Nitride-Oxide-Coated versus Everolimus-Eluting Stents in acute coronary syndrome the randomized TIDES-ACS-Trial, JACC Cardiovascular Intervention
- (2) López-Minguez JR, et al, TIOMAX; A Spanish Multicenter Registry of the real-world use of the Titanium Optimax biostent, Catheter Cardiovascular Intervention
- (3) Kajalainen PP et al, Clinical outcome of titanium-nitride-oxide-coated cobalt chromium stents in patients with de novo coronary lesions: 12-month results of the OPTIMAX first-in-man study, Catheter Cardiovascular Intervention

5.2 Potential Adverse Effects

Potential adverse effects associated with the stenting procedure include but are not limited to:

- Allergic reaction to contrast media, antiplatelet aggregation and/or anticoagulant medications, stent TiNO coating
- Angina
- Arteriovenous fistula
- Arrhythmias, including ventricular fibrillation
- Bleeding
- Coronary artery spasm
- Death
- Embolism
- Fever
- Haemorrhage or haematoma
- Hypotension/hypertension
- Infection, local or systemic (sepsis)
- Myocardial infarction (MI)
- Need for immediate coronary bypass surgery
- Vessel and/or coronary artery injury, perforation, rupture and/or dissection
- Restenosis of the dilated artery
- Stent migration
- Stent thrombosis
- Vessel and/or coronary artery total occlusion

-VI- WARNINGS

- The TITAN OPTIMAX stent should only be implanted in hospitals or clinics authorised to perform coronary angioplasty and by physicians with technical expertise in coronary angioplasty who may be assisted by physicians capable of treating major clinical complications with cardiac surgery.
- Check the expiry date on the packaging protecting the product's sterility. Do not use the device if the expiry date has passed.

- Do not re-sterilise, reprocess and/or re-use the device.
- Do not re-use the stent if it fails to pass. Re-sterilisation, reprocessing and/or re-use may compromise the performance of the device and its integrity. Such actions may result in contamination of the device and/or cause patient infection or cross-infection. HEXACATH cannot be held responsible for any accidental, direct or consequential damage resulting from re-sterilisation or re-use of the product.
- Ensure the integrity of the catheter before use.
- Do not use if the stent or catheter is damaged (leakage, breakage, cracking, loose stent, stent not centred on the balloon etc.).
- Observe aseptic conditions during all phases of use.
- Do not use gauze pads; the fibres can damage the stent.
- Do not loosen and re-tighten the stent on another catheter.
- Do not reposition the stent on its catheter.
- Do not use any gas, air or other liquids to inflate the balloon other than those recommended in the instruction manual.
- Never inflate the balloon before final stent placement at a pressure exceeding 0.5 bar. This could result in the stent opening prematurely and being unable to progress through the artery.
- Do not inflate the balloon beyond the rated burst pressure (RBP). This increases the risk of balloon rupture, which can lead to vessel occlusion, balloon entrapment and associated complications.
- The catheter must be observed under fluoroscopy throughout the procedure.
- Never move the balloon with a pre-mounted stent forwards without ensuring that the stent is perfectly attached to the balloon. This could result in the loss of the stent before it is deployed in the vessel and could lead to acute myocardial infarction or death.
- Stents can cause nuclear magnetic resonance (or MRI) artefacts due to magnetic field distortion. These artefacts caused by the cobalt-chrome alloy of the stent are comparable to those caused by metallic surgical clips. In order to minimise the risk of stent migration under a strong magnetic field, it is recommended that a nuclear magnetic resonance examination is performed only after complete endothermalization of the stent, i.e. at least three months after implantation.

-VII- PRECAUTIONS

- Use the device before the expiry date.
- Store above 0°C and below 40°C, protected from light and moisture.
- Do not use if package has been damaged or opened.
- Before performing the angioplasty, check that the device is working properly and that it is the right size and shape for the procedure. Do not use the device if there is any damage or if there is any doubt about its integrity.
- Prior to inserting the catheter, administer the appropriate dose of coronary anticoagulant and vasodilator.
- Check that the air exhaust in each system is complete and that there are no leaks in the various connections.

-VIII- RECOMMENDED EQUIPMENT

- A puncture kit
- Appropriately sized guide catheter: 5F guide catheters, with a minimum inner lumen of 0.058" (1.47 mm), allow for the advancement of pre-mounted TITAN OPTIMAX stents from 2.0 mm to 4.0 mm in diameter. The 6F guide catheters, with a minimum internal lumen of 0.064" (1.62 mm), allow for the advancement of pre-mounted TITAN OPTIMAX stents from 4.5 and 5.0 mm in diameter).
- A 0.014" coronary guide of appropriate length
- An appropriate syringe
- An appropriately sized haemostasis valve
- A 3-way stopcock
- An inflation device
- Heparinized saline solution
- Contrast agent diluted with saline solution in a 1:1 ratio
- A flushing needle

-IX- OPERATING INSTRUCTIONS**9.1 Preparation of the vessel to be treated**

Conventional angioplasty techniques for placement of an arterial introducer, guide catheter and guidewire should be used. The guide catheter must have a radiopaque distal end. The haemostasis valve must be able to accept a 5F or 6F guide catheter. The recommended guidewire should be 0.014".

9.2 Stent size selection

Longer stents require a larger radius of curvature and a larger vessel diameter to pass angled sites. The length of the stent must correspond at least to the length of the lesion to be treated. The use of small stents is recommended over the use of long stents if access to the lesion is relatively narrow or tortuous or if there is any doubt about ease of access. Stent diameter selection: a suitable stent diameter should be selected, based on the recommended stent diameter to artery diameter ratio of 1.1:1.

9.3 Preparation of the installation system

1. Make a solution of 50% contrast agent and 50% saline solution.
 2. Carefully remove the installation system (catheter) from its plastic packaging. Remove the protective tube from the distal end of the system, exposing the stent to the balloon. Check that the folded stent has a uniform surface, is free of protrusions, and is properly centred on the balloon (between the two radiopaque markers). In the event of a problem, return the product to HEXACATH.
 3. Connect a suitable syringe filled with approximately 10ml of saline solution to an appropriate flushing needle. Insert the needle into the distal end of the catheter and flush the guidewire lumen with low pressure. Remove the syringe and needle.
- Avoid manipulating the stent while flushing the guide lumen, as this may alter the position of the stent on the balloon. Be especially careful not to change the stent's position on the balloon. This is particularly important when removing the delivery system from its packaging, placing it on the guidewire and inserting it into the Y-haemostasis valve (internal diameter > 1.67mm) and the hub of the guide catheter.

4. Connect a 3-way stopcock to the catheter luer-lock port.
 5. Fill the syringe with 3 ml of the solution of 50% contrast agent 50% saline solution. Drain any remaining air and attach it to the 3-way stopcock.
 6. Open the 3-way stopcock.
 7. Hold the syringe and the proximal end of the installation system in an upright position, pointing the plunger of the syringe upwards.
 8. Apply negative pressure by pulling the plunger of the syringe out as far as possible. Maintain the negative pressure until all air bubbles have disappeared from the solution.
- If negative pressure cannot be maintained, check the connections. If the problem persists, return the product to HEXACATH.
9. Turn the 3-way stopcock to the closed position at the inflation port, then remove the syringe and purge any remaining air from the syringe.
 10. Reattach the syringe to the inflation port, open the 3-way stopcock and repeat steps 4 and 5.
 11. Before removing the syringe, release the plunger and allow the contrast agent to flow into the inflation lumen and balloon.
 12. Close the 3-way stopcock at the inflation port, then remove the syringe.
 13. Fill the inflation device with around 3 ml of the solution of 50% contrast agent 50% saline solution. Purge any remaining air from the system.
 14. Connect the inflation device without introducing air and attach it to the balloon port using the 3-way stopcock.
 15. Immerse the distal portion of the catheter with the hydrophilic coating in sterile saline solution before use.
 16. Insert the introducer and then the appropriately sized guide catheter.
 17. Under fluoroscopy, insert and advance the coronary guidewire through the guide catheter until the desired location is reached.

9.4 Insertion of the TITAN OPTIMAX stent

1. Insert the TITAN OPTIMAX stent and installation system into the Y-hemostatic valve (I.D. > 1.67 mm) by inserting the proximal part of the guide catheter (0.014" max.) at the catheter tip. Ensure that the valve is fully open when the stent is inserted. Ensure that the stent is placed in the center of the balloon of the delivery stent system. If not, remove the entire system and repeat the procedure with a new device. If resistance is encountered, do not force the passage, as the stent may be damaged or displaced on the balloon.
2. Advance the stent delivery system into the guide catheter under fluoroscopic guidance. The stent can be used to determine when the end of the balloon reaches the end of the guide catheter. The output markers can also be used to determine the system's position. Follow recommended PTCA techniques to move the delivery system towards the lesion

Avoid over-tightening the valve, as this may obstruct the flow of contrast agent to the balloon, thus slowing down inflation and deflation.

3. While observing the TITAN OPTIMAX radiopaque stent by fluoroscopy, place the stent in the lesion. Before deploying the stent, perform a high-definition fluoroscopic examination to ensure that the stent has not been damaged or displaced during positioning. If the initial position of the guide catheter is lost, avoid pulling or pushing it over the stent. Movement of the distal end of the guide catheter may damage the stent.
4. Perform one or more inflations to achieve a diameter relative to the artery of 1.1/1. Check the stent expansion before removing the balloon to minimise shrinkage. Evaluate this ratio once the stent is fully deployed, using orthogonal angiographic visualisation. The angiographic image of the stent should appear larger than the vessel itself. The measured diameter of the stent should be approximately 10% larger than the vessel diameter before stent insertion. Intravascular ultrasound may be useful for assessing stent expansion. Using a stent that is too large and applying a higher inflation pressure than recommended can cause distal dissection. The nominal burst pressure (RBP) must not be exceeded under any circumstances. Do not deploy the stent if it is not properly placed in the stenotic segment of the vessel. If its position is not optimal, the stent must be moved or removed. Balloon inflation should be performed with the guidewire extending beyond the end of the catheter. It is strongly recommended that the guidewire, balloon catheter, or both be left in place beyond the lesion until the procedure is complete.

9.5 General recommendations

The stent and in particular its ends must be perfectly adherent to the balloon otherwise its advancement could be compromised.

9.6 Specific recommendations**9.6.1. Placement of multiple stents**

The following recommendations should be considered if the lesion(s) requires more than one stent. The first stent should be placed distal to the lesion and subsequent stents should be placed from farthest to closest. Commonly in lesions that require multiple stents, they must be placed side by side or slightly overlapping to stabilize the dissections or to avoid wall protrusions between the stents. If an initial angiographic result is imperfect the stent should be progressively oversized with larger diameter balloon catheters.

9.6.2. Removal and disassembly

After opening the stent, deflate the balloon, allowing sufficient time for complete deflation. Very slowly withdraw the balloon from the stent, maintaining negative pressure. Hold the guide catheter in place to prevent over-intubation of the guide catheter. The balloon must be completely deflated and checked by fluoroscopic inspection before removal. Remove the deflated delivery system, while keeping the guidewire and guide catheter in position in the vessel.

9.6.3. Removal of an unopened stent

During access to the lesion, if the progression of the stent encounters unusual resistance, it is recommended that the guide catheter and the installation system are removed in bloc. Excessive thrust on the delivery system may result in damage or loss of the stent. In some cases, the stent and its catheter may not be advanced towards the target lesion, in which case they must be removed. Removal of a stent from an undeployed balloon may result in the loss of the stent in the coronary artery. In the event of unusual resistance, the stent and balloon should be removed together with the guide catheter and guidewire to avoid loss of the TITAN OPTIMAX stent in the coronary system. The guide catheter, stent and balloon should ideally be removed in one piece through the introducer under fluoroscopic guidance. However, if it is felt that there is a risk of losing the stent through the introducer (deformed or compacted stent), it is possible to ensure passage by using a recovery lasso to fix the stent, and/or by exchanging the 5F/6F introducer for an 8F.

9.7 Magnetic resonance imaging (MRI)

Non-clinical testing and MRI simulations were performed to evaluate the entire family of the TITAN OPTIMAX Coronary Artery Stent. Non-clinical testing demonstrated that the entire family of this product TITAN OPTIMAX Coronary Artery Stent is MR Conditional. A patient with an implant from this family can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla or 3-Tesla, only.
- Maximum spatial gradient magnetic field of 1,900-gauss/cm (19-T/m).

- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode. Under the scan conditions defined, an implant from the TITAN OPTIMAX Coronary Artery Stent is expected to produce a maximum temperature rise of 3.1°C after 15-minutes of continuous scanning (i.e., per pulse sequence). In non-clinical testing, the image artifact caused by and implant from the TITAN OPTIMAX Coronary Artery Stent extends approximately 5-mm from this device when imaged with a gradient echo pulse sequence and a 3-Tesla MR system.

9.8 Disposal procedure

After use, dispose of the product and its packaging in accordance with local regulations to avoid any risk of infection by substances of human origin.

X- CLINICAL BENEFIT, SAFETY AND PERFORMANCE**10.1 Performance feature**

TITAN OPTIMAX coronary stent has the ability to reopen stenotic arteries when deployed and minimise post-operative restenosis (due to its specific TiNO coating).

10.2 Clinical benefit

The main clinical benefits expected from performance of TITAN OPTIMAX coronary stent is to prevent a heart attack (myocardial infarction) or to minimise the consequences of a heart attack. As well, stenting is a minimally invasive procedure which does not require general anaesthesia and major surgery for the patient.

10.3 Summary of safety and performance features

The summary of safety and performance features can be found at: www.hexacath.com

-XI- PRODUCT LIMITATIONS AND WARRANTIES

HEXACATH guarantees that the TITAN OPTIMAX system has been designed, manufactured, tested and packaged with the utmost care. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. The handling and storage of this device, as well as factors relating to the patient, diagnosis, treatment, surgical procedures and other areas outside of HEXACATH's control, directly affect the device and the results obtained from its use. HEXACATH's obligations are limited to device replacement. HEXACATH shall not be liable in any way for any loss, damage or expense, incidental or consequential, arising from the use of the device. HEXACATH does not assume, nor does it authorise any third party to assume on its behalf, any additional liability or obligation in connection with the device. HEXACATH shall not be liable for any re-use, reprocessing or re-sterilisation of the products and makes no warranties, express or implied, including but not limited to any warranty of merchantability or fitness for a particular purpose with respect to these products.

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- (1) Tonino P. et al., Titanium-Nitride-Oxide-Coated versus Everolimus-Eluting Stents in acute coronary syndrome the randomized TIDES-ACS-Trial, JACC Cardiovascular Intervention
- (2) López-Minguez JR. et al, TIOMAX, A Spanish Multicenter Registry of the real-world use of the Titanium Optimax bioactive, Catheter Cardiovascular Intervention
- (3) Karjalainen PP. et al, Clinical outcome of titanium-nitride-oxide-coated cobalt chromium stents in patients with de novo coronary lesions: 12-month results of the OPTIMAX first-in-man study, Catheter Cardiovascular Intervention

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- (1) Tonino P. et al., Titanium-Nitride-Oxide-Coated versus Everolimus-Eluting Stents in acute coronary syndrome the randomized TIDES-ACS-Trial, JACC Cardiovascular Intervention
- (2) López-Minguez JR. et al, TIOMAX; A Spanish Multicenter Registry of the real-world use of the Titanium Optimax biostent, Catheter Cardiovascular Intervention
- (3) Karjalainen PP. et al, Clinical outcome of titanium-nitride-oxide-coated cobalt chromium stents in patients with de novo coronary lesions: 12-month results of the OPTIMAX first-in-man study, Catheter Cardiovascular Intervention

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- (1) Tonino P, et al., *Titanium-Nitride-Oxide-Coated versus Everolimus-Eluting Stents in acute coronary syndrome the randomized TIDES-ACS-Trial*, JACC Cardiovascular Intervention
- (2) López-Minguez JR, et al., *TIOMAX; A Spanish Multicenter Registry of the real-world use of the Titanium Optimax bioestent, Catheter Cardiovascular Intervention*
- (3) Karjalainen PP, et al., *Clinical outcome of titanium-nitride-oxide-coated cobalt chromium stents in patients with de novo coronary lesions: 12-month results of the OPTIMAX first-in-man study, Catheter Cardiovascular Intervention*

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- (1) Tonino P. et al., Titanium-Nitride-Oxide-Coated versus Everolimus-Eluting Stents in acute coronary syndrome the randomized TIDES-ACS-Trial, JACC Cardiovascular Intervention
- (2) López-Minguez JR. et al, TiOMAX, A Spanish Multicenter Registry of the real-world use of the Titanium Optimax bioestent, Catheter Cardiovascular Intervention
- (3) Karjalainen PP. et al, Clinical outcome of titanium-nitride-oxide-coated cobalt chromium stents in patients with de novo coronary lesions: 12-month results of the OPTIMAX first-in-man study, Catheter Cardiovascular Intervention

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- (1) Tonino P. et al., Titanium-Nitride-Oxide-Coated versus Everolimus-Eluting Stents in acute coronary syndrome the randomized TIDES-ACS-Trial, JACC Cardiovascular Intervention
- (2) López-Minguez JR. et al, TIOMAX: A Spanish Multicenter Registry of the real-world use of the Titanium Optimax biostent, Catheter Cardiovascular Intervention
- (3) Karjalainen PP. et al, Clinical outcome of titanium-nitride-oxide-coated cobalt chromium stents in patients with de novo coronary lesions: 12-month results of the OPTIMAX first-in-man study, Catheter Cardiovascular Intervention

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- (1) Tonino P. et al., *Titanium-Nitride-Oxide-Coated versus Everolimus-Eluting Stents in acute coronary syndrome the randomized TIDES-ACS-Trial*, JACC Cardiovascular Intervention
- (2) López-Minguez JR. et al., *TIOMAX; A Spanish Multicenter Registry of the real-world use of the Titanium Optimax biostent, Catheter Cardiovascular Intervention*
- (3) Karjalainen PP. et al., *Clinical outcome of titanium-nitride-oxide-coated cobalt chromium stents in patients with de novo coronary lesions: 12-month results of the OPTIMAX first-in-man study*, Catheter Cardiovascular Intervention

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- (1) Tonino P. et al., Titanium-Nitride-Oxide-Coated versus Everolimus-Eluting Stents in acute coronary syndrome the randomized TIDES-ACS-Trial, JACC Cardiovascular Intervention
- (2) López-Minguez JR. et al., TIOMAX; A Spanish Multicenter Registry of the real-world use of the Titanium Optimax bioent, Catheter Cardiovascular Intervention
- (3) Karjalainen PP. et al., Clinical outcome of titanium-nitride-oxide-coated cobalt chromium stents in patients with de novo coronary lesions: 12-month results of the OPTIMAX first-in-man study, Catheter Cardiovascular Intervention

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(1) Tonino P. et al., Titanium-Nitride-Oxide-Coated versus Everolimus-Eluting Stents in acute coronary syndrome the randomized TIDES-ACS-Trial, JACC Cardiovascular Intervention

(2) López-Minguez JR. et al., TIOMAX: A Spanish Multicenter Registry of the real-world use of the Titanium Optimax bioactive, Catheter Cardiovascular Intervention

(3) Karjalainen PP. et al., Clinical outcome of titanium-nitride-oxide-coated cobalt chromium stents in patients with de novo coronary lesions: 12-month results of the OPTIMAX first-in-man study, Catheter Cardiovascular Intervention

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- (1) Tonino P. et al., *Titanium-Nitride-Oxide-Coated versus Everolimus-Eluting Stents in acute coronary syndrome the randomized TIDES-ACS-Trial, JACC Cardiovascular Intervention*
- (2) López-Minguez JR, et al., *TIOMAX, A Spanish Multicenter Registry of the real-world use of the Titanium Optimax biostent, Catheter Cardiovascular Intervention*
- (3) Karjalainen PP, et al., *Clinical outcome of titanium-nitride-oxide-coated cobalt chromium stents in patients with de novo coronary lesions: 12-month results of the OPTIMAX first-in-man study, Catheter Cardiovascular Intervention*

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- (1) Tonino P. et al., Titanium-Nitride-Oxide-Coated versus Everolimus-Eluting Stents in acute coronary syndrome the randomized TIDES-ACS-Trial, JACC Cardiovascular Intervention
- (2) López-Minguez JR. et al, TIOMAX; A Spanish Multicenter Registry of the real-world use of the Titanium Optimax bioactive, Catheter Cardiovascular Intervention
- (3) Kajalainen PP. et al, Clinical outcome of titanium-nitride-oxide-coated cobalt chromium stents in patients with de novo coronary lesions: 12-month results of the OPTIMAX first-in-man study, Catheter Cardiovascular Intervention

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(1) Tonino P et al., *Titanium-Nitride-Oxide-Coated versus Everolimus Eluting Stents in acute coronary syndrome the randomized TIDES-ACS-Trial, JACC Cardiovascular Intervention*

(2) López-Minguez JR, et al., *TIOMAX: A Spanish Multicenter Registry of the real-world use of the Titanium Optimax biostent, Catheter Cardiovascular Intervention*

(3) Karjalainen PP, et al., *Clinical outcome of titanium-nitride-oxide-coated cobalt chromium stents in patients with de novo coronary lesions: 12-month results of the OPTIMAX first-in-man study, Catheter Cardiovascular Intervention*

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- (1) Tonino P, et al., Titanium-Nitride-Oxide-Coated versus Everolimus Eluting Stents in acute coronary syndrome the randomized TIDES-ACS-Trial, JACC Cardiovascular Intervention
- (2) López-Minguez JR, et al, TIOMAX; A Spanish Multicenter Registry of the real-world use of the Titanium Optimax bioactive, Catheter Cardiovascular Intervention
- (3) Kajalainen PP, et al, Clinical outcome of titanium-nitride-oxide-coated cobalt chromium stents in patients with de novo coronary lesions: 12-month results of the OPTIMAX first-in-man study, Catheter Cardiovascular Intervention

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- (1) Tonino P. et al., Titanium-Nitride-Oxide-Coated versus Everolimus-Eluting Stents in acute coronary syndrome the randomized TIDES-ACS-Trial, JACC Cardiovascular Intervention
- (2) López-Minguez JR. et al, TiOMAX, A Spanish Multicenter Registry of the real-world use of the Titanium Optimax bioinstent, Catheter Cardiovascular Intervention
- (3) Karjalainen PP. et al, Clinical outcome of titanium-nitride-oxide-coated cobalt chromium stents in patients with de novo coronary lesions: 12-month results of the OPTIMAX first-in-man study, Catheter Cardiovascular Intervention

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