Mistral Pro

SC

Mistral Pro

NC

HEXACATH

PIONEER IN BIO ACTIVE COATING

CATHÉTER DE DILATATION POUR APTC MISTRAL PRO SC / MISTRAL PRO NC	
NOTICE D'UTILISATION - Français	2
DILATATION CATHETER FOR APTC MISTRAL PRO SC / MISTRAL PRO NC	
INSTRUCTIONS FOR USE - English	4
DILATATIONS KATHETER FÜR APTC MISTRAL PRO SC / MISTRAL PRO NC	
GEBRAUCHSANWEISUNG - Deutsch	6
CATETERE DI DILATAZIONE PER APTC MISTRAL PRO SC / MISTRAL PRO NC	
MANUALE UTENTE - Italiano	9
CATÉTER DE DILATACIÓN PARA APTC MISTRAL PRO SC / MISTRAL PRO NC	
MANUAL DEL USUARIO - Español	11
CATETER DE DILATAÇÃO PARA APTC MISTRAL PRO SC / MISTRAL PRO NC	
MANUAL DE UTILIZAÇÃO - Português	13
PTCA 冠脉球囊扩张导管 MISTRAL PRO SC / MISTRAL PRO NC	
使用说明书 - 中國	15
ΚΑΘΗΤΗΡΑΣ ΔΙΑΣΤΟΛΗΣ MISTRAL PRO SC / MISTRAL PRO NC ΓΙΑ PTCA	
ΕΓΧΕΙΡΙΔΙΟ ΧΡΗΣΤΗ - Ελληνικά	17
CEWNIK ROZSZERZAJĄCY DO PTCA MISTRAL PRO SC / MISTRAL PRO NC	
INSTRUKCJA UŻYTKOWANIA - Polski	20
DILATAČNÍ KATÉTR PRO APTC MISTRAL PRO SC / MISTRAL PRO NC	
UŽIVATELSKÁ PŘÍRUČKA - Český	22
DILATAČNÝ KATÉTER PRE PTCA MISTRAL PRO SC/MISTRAL PRO NC	
INFORMÁCIE PRE POUŽÍVATEĽOV - Slovenský	24
DILATATIONS KATETER TIL APTC MISTRAL PRO SC / MISTRAL PRO NC	
BRUGERMANUAL - Dansk	26
DILATASJONSKATETER FOR APTC MISTRAL PRO SC / MISTRAL PRO NC	
BRUKSINSTRUKSJONER - Norsk	28
DILATATIONS KATETER FÖR APTC MISTRAL PRO SC / MISTRAL PRO NC	
BRUKSANVISNING - Svenska	30
DILATÁCIÓS KATÉTER AZ APTC MISTRAL PRO SC / MISTRAL PRO NC ESZKÖZÖKHÖZ	
HASZNÁLATI ÚTMUTATÓ - Magyar	32
CATHETER NONG BÓNG CHO PTCA MISTRAL PRO SC / MISTRAL PRO NC	
HƯỚNG DẪN SỬ DỤNG - Tiếng Việt	34

Recyclable paper

REV04-06/2023 - PEA-06/2023

Marquage
CE
2022

CE
0459

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DILATATION CATHETER FOR APTC MISTRAL PRO SC / MISTRAL PRO NC INSTRUCTIONS FOR USE

Read all instructions carefully before use.

Follow all warnings and precautions in the instructions to avoid any possible complications.

Any serious incident occurring in connection with the device must be notified to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

Sterile. Sterilised with Ethylene Oxide. Non-pyrogenic. Single use. Do not reuse, reprocess or resterilise. Do not use if the packaging has been damaged or opened. Use before the expiry date.

After use, dispose of the product and packaging in accordance with local regulations.

English

DILATATION CATHETER FOR APTC MISTRAL PRO SC / MISTRAL PRO NC manual
©2022 HEXACATH
All rights reserved
MISTRAL PRO SC / PRO NC, Hexacath and their respective logos are registered trademarks of Hexacath S.A.S

-I- DESCRIPTION

The MISTRAL PRO SC and MISTRAL PRO NC catheters are intended to be used in adult patient with clinical symptoms of myocardial ischemia related to the pathological condition of coronary arteries in order to re-open the artery to its native diameter at the level of a lesion, following angioplasty.

The MISTRAL PRO SC and MISTRAL PRO NC catheters are Percutaneous Transluminal Coronary Angioplasty (PTCA) catheters with a rapid exchange design and a balloon at the tip.

The material used gives the balloon a controlled compliance allowing the diameter to be known in accordance with the inflation pressure (MISTRAL PRO SC: semi-compliant balloon, MISTRAL PRO NC: non-compliant balloon). Depending on its length, the balloon has one or two radiopaque markers for correct positioning under fluoroscopy.

The distal part of the catheter (25cm), coated with a hydrophilic coating, has a double coaxial lumen. The outer lumen is used to inflate the balloon and the inner lumen allows the introduction of a coronary guidewire to facilitate the catheter's progression towards the desired dilatation site. The catheter is fitted with a tapered tip to facilitate its progress towards and through the stenosis. The proximal part of the catheter consists of a stainless steel hypotube with a luer lock connector at the tip for inflating and deflating the balloon. Two markers located 90 cm (brachial/radial approach) and 100 cm (femoral approach) from the distal end of the catheter allow it to be positioned at the end of the guide catheter without the need for fluoroscopy.

A protective tube is positioned over the balloon and a mandrel is inserted into the inner lumen to protect the integrity of the device.

-II- PRESENTATION AND REFERENCES

Packaged in a peelable pouch (sterile barrier system) and cardboard box.

Contents: One (1) MISTRAL PRO NC or MISTRAL PRO SC catheter, one (1) instruction leaflet and one (1) compliance card.

The MISTRAL PRO SC catheter is available in lengths of 10, 15, 20, 30, 40 mm and diameters of 1.25, 1.5, 2.0, 2.25, 2.5, 2.75, 3.0, 3.5, 4.0, 4.5, 5.0 mm.

The MISTRAL PRO NC catheter is available in lengths of 8, 11, 15, 20, 30 mm and diameters of 2.0, 2.25, 2.5, 2.75, 3.0, 3.25, 3.5, 3.75, 4.0, 4.5, 5.0 mm.

-III- INDICATIONS

MISTRAL PRO SC and MISTRAL PRO NC catheters should only be used by physicians trained in PTCA procedures. They are designed for use in adult patients with manifestations of myocardial ischaemia related to coronary artery disease and who are candidates for myocardial revascularisation.

MISTRAL PRO SC and MISTRAL PRO NC catheters are indicated for:

- Balloon dilatation of the stenotic portion of a coronary artery stenosis, for the purpose of improving myocardial perfusion

- Balloon dilatation of a stent after implantation.

The MISTRAL PRO NC catheter should preferably be used for expansion after stenting (post-dilatation), as well as for calcified and/or fibrous lesions.

-IV- CONTRAINDICATIONS

- Patients with chronic total occlusion of the target coronary artery
- Coronary artery spasm in the absence of a significant stenosis
- Severe stenosis of the unprotected left main coronary artery
- Allergy to contrast media

-V- ADVERSE EFFECTS

5.1 Observed adverse effects

Observed adverse effects comes from a systematic review of the literature. They are summarized in the table below:

Adverse events	Quantification (%) – [Occurrence/Total number of subjects involved in the study]
Haemorrhage / haematoma	1,96% ^{1,2} - [1/51]
Unstable angina	1,7% ³ (non-compliant balloons) - [1/59] 3,3% ⁴ (compliant balloons) - [2/60]
Acute myocardial infarction	1,96% ^{1,2} - [1/51] 6,7% ^{5,5} - [4/60]
Pseudoaneurysm	1,67% ^{5,5} - [1/60]

¹ *"Sprinter Legend Balloon Catheter - Study ResultsClinicalTrials.Gov"*, <https://clinicaltrials.gov/ct2/show/results/NCT00961311>.

² David E. Kandzari et al., "Clinical Outcomes Following Predilatation with a Novel 1.25-Mm Diameter Angioplasty Catheter", *Catheterization and Cardiovascular Interventions: Official Journal of the Society for Cardiac Angiography & Interventions* 77, n° 4 (1 mars 2011): 510-14, <https://doi.org/10.1002/ccd.22734>

³ Erdem Özel et al., "What Is Better for Predilatation in Bioresorbable Vascular Scaffold Implantation: A Non-Compliant or a Compliant Balloon?", *Anatolian Journal of Cardiology* 16, n° 4 (avril 2016): 244-49, <https://doi.org/10.5152/AnatolJCardiol.2015.6184>

⁴ https://www.bostonscientific.com/content/dam/Manuals/us/current-rev-en/90990229-01A_Emerge_dtu_en-US_S.pdf

⁵ "Evaluation of Coronary Luminal Diameter Enlargement With Emerge™ 1,20 mm PTCA Dilatation Catheter - Study Results - ClinicalTrials.Gov", <https://clinicaltrials.gov/ct2/show/results/NCT01635881>

5.2 Potential Adverse Effects

The following adverse events may be a consequence of a coronary angioplasty procedure:

- Infection
- Haemorrhage / haematoma
- Embolism
- Thrombosis
- Coronary spasms
- Coronary vessel injury (dissection, perforation, rupture)
- Arteriovenous fistula
- Restenosis of the dilated artery
- Total occlusion
- Low/high blood pressure

- Unstable angina
- Arrhythmias, including ventricular fibrillation
- Acute myocardial infarction
- Transient ischaemic attacks
- Death
- Need for immediate coronary bypass surgery
- Pseudoaneurysm

-VI- WARNINGS

- These devices are sterilised with ethylene oxide and are intended for single use. Do not reuse, reprocess or resterilise. Reprocessing, re-sterilisation and/or re-use of these devices may compromise their performance and integrity. Such actions may result in contamination of the device and/or cause patient infection or cross infection.
- MISTRAL PRO SC and MISTRAL PRO NC catheters should only be used in hospitals or clinics licensed to perform coronary angioplasty and in facilities that treat serious or life-threatening complications. Observe aseptic conditions during all phases of use.
- Inside the vascular system, the catheter must be observed under fluoroscopy throughout the procedure. The catheter should only be pushed or removed when the balloon is fully deflated. If any resistance is felt, the cause should be determined before proceeding. Moving the catheter in an attempt to overcome the resistance may result in vascular injury and damage to the catheter.
- Do not use or attempt to straighten a catheter, the body of which has been bent or twisted. This could cause the catheter body to rupture.
- The inflation pressure of the device must not exceed the balloon's rated burst pressure (RBP) shown on the compliance card. Excessive inflation pressure can cause the balloon to rupture.
- The diameter of the inflated balloon should not exceed the diameter of the vessel proximal and distal to the stenosis, to minimise the risk of vascular damage.
- Do not use air or gas to inflate the balloon.
- PTCA in patients who are not acceptable candidates for coronary artery bypass graft surgery requires careful consideration, including possible hemodynamic support during PTCA as treatment of this patient population carries special risk.
- Use extreme caution and careful judgment in patients for whom anticoagulation is not indicated and in patients with impaired renal function.

-VII- PRECAUTIONS

- Use before the expiry date.
- Store above 0°C and below 40°C, away from light and moisture.
- Do not use if the packaging has been damaged or opened.
- Before performing the angioplasty, check that the device is working properly and that it is the right size and shape for the procedure. Do not use the device if there is any damage or if there is any doubt about its integrity.
- Before inserting the catheter, administer the appropriate dose of coronary anticoagulant and vasodilator.

- Check that the air venting in each system is complete and that there are no leaks in the various connections.
- Resistance may be felt with a MISTRAL PRO during insertion into or withdrawal from the guide catheter. If this is the case, choose a larger guide catheter. For MISTRAL PRO with a diameter of 4.0 mm and a length of 40 mm and for MISTRAL PRO with a diameter of 4.5 and 5.0 mm, use a guide catheter with a minimum diameter of 6F.

-VIII- RECOMMENDED EQUIPMENT

- A puncture kit
- A guide catheter with a minimum diameter of 5F (0.058")
- A guide catheter with a minimum diameter of 6F (0.071") for MISTRAL PRO with a diameter of 4.0 mm and a length of 40 mm and for MISTRAL PRO with a diameter of 4.5 and 5.0 mm
- A 0.014" coronary guidewire of appropriate length
- A suitable syringe
- An appropriately sized haemostasis valve
- A 3-way valve
- An inflation device (recommended with a pressure gauge)
- Contrast medium diluted with saline in a 1:1 ratio
- A flushing needle

-IX- INSTRUCTIONS FOR USE

9.1. Preparation

Preparing the inflation device

- Prepare the inflation device according to the manufacturer's instructions.
- Purge the air from the system.

Preparing the MISTRAL PRO NC or SC catheter

- Make a solution of 50% contrast medium and 50% saline.
- Carefully remove the MISTRAL PRO catheter from its packaging and the spiral.
- Keep the shield on the balloon and the mandrel in the inner lumen of the catheter.
- Connect a 3-way valve to the MISTRAL PRO catheter hub.
- Draw up 1-3 ml of 50/50 contrast medium solution into a suitable syringe.
- Connect the filled syringe to the 3-way valve and open the valve.
- Hold the tip of the syringe down and create a strong vacuum for 15 to 20 seconds until no more bubbles appear. Slowly release the piston. Close the tap and remove the syringe. Evacuate all the air contained in the syringe.
- Repeat the operation if necessary.
- Connect the inflation device without introducing air into the system and open the valve. Apply negative pressure and close the valve.
- Connect a suitable syringe filled with approximately 10 ml of saline to a suitable flushing needle.
- Remove the balloon protection and the mandrel. If there is any unusual resistance when removing the balloon protection or the mandrel, the MISTRAL PRO catheter should not be used and should be replaced by another one.
- Insert the needle into the distal end of the MISTRAL PRO catheter and bleed the lumen of the guidewire with low pressure.
- Immerse the distal part of the MISTRAL PRO catheter covered with the hydrophilic coating (for about 25 cm) in sterile saline before use.

9.2. Procedure for use

- Connect a haemostasis valve to the appropriate guide catheter hub that has been placed in the vascular system.
- Insert the coronary guidewire into the guide catheter through the haemostasis valve.

- Under fluoroscopy, advance the coronary guidewire to the relevant vessel and place it in the desired location.
- Insert the proximal portion of the coronary guidewire into the distal end of the MISTRAL PRO catheter and monitor its exit through the proximal port.
- Loosen the haemostasis valve and advance the MISTRAL PRO catheter to the distal end of the guidewire while holding the coronary guidewire stationary. The two markers on the hypotube allow, depending on a brachial, radial or femoral approach, the position of the MISTRAL PRO catheter in the guide catheter to be identified.
- Under fluoroscopy, use radiopaque markers for balloon positioning in the stenosis or in the stent for MISTRAL PRO NC. If the stenosis cannot be crossed with the selected MISTRAL PRO catheter, use a smaller diameter catheter to pre-dilate the stenosis and facilitate the passage of a more appropriate sized catheter. Do not move the MISTRAL PRO catheter without a coronary guidewire in front of it.
- Dilate the stenosis or stent for MISTRAL PRO NC using a recognised coronary angioplasty technique to a pressure below the balloon's rated burst pressure (RBP). If there is any difficulty in inflating the balloon, do not continue with the MISTRAL PRO catheter and remove it. In case of post-dilatation with the MISTRAL PRO NC, it is recommended to select a balloon shorter than the stent in order to avoid trauma to the arterial wall outside the stent.
- Perform balloon dilation with the coronary guidewire extending beyond the distal end of the MISTRAL PRO catheter. It is recommended that the coronary guidewire and/or balloon be left in place through the stenosis or through the stent for MISTRAL PRO NC until the procedure is complete.
- After angioplasty or post-stent dilation, deflate the balloon and remove the MISTRAL PRO catheter from the lesion.
- Maintain the coronary guidewire in the dilated artery and ensure the result using a recognised coronary angiography technique.
- After the angiogram, slowly remove the coronary guidewire and deflated balloon from the guide catheter.

9.3. Exchange procedure

- Loosen the haemostasis valve, hold the coronary guidewire and valve in one hand while the other hand holds the body of the MISTRAL PRO SC or PRO NC catheter.
- Hold the coronary guidewire in place and begin pulling on the deflated catheter, monitoring the position of the coronary guidewire under fluoroscopy during the exchange.
- Continue to pull until the proximal access valve of the coronary guidewire tract appears, then proceed as for a coaxial catheter.
- Prepare the new catheter as described in the section Preparing the catheter.
- Insert the new catheter as described in the Procedure for use section.

9.4 Disposal procedure

After use, dispose of the product and its packaging in accordance with local regulations to avoid the risk of infection with substances of human origin.

-X- CLINICAL BENEFITS, SAFETY AND PERFORMANCE

10.1. Performance characteristics

- Ability to access the lesion
- Ability to open the stenotic artery
- Ability to optimise stent deployment during post-dilatation (MISTRAL PRO NC)

10.2. Clinical benefits

The expected clinical benefit of the use of a MISTRAL PRO catheter is to help the restoration of arterial circulation, permitting to prevent a heart attack (myocardial infarction) or to minimize the consequences of heart attack. As well, balloon angioplasty is a minimally invasive procedure which does not require general anesthesia and major surgery.

10.3 Summary of safety and performance characteristics

The summary of safety and performance characteristics can be found at: www.hexacath.com

-XI- WARRANTY

HEXACATH guarantees that the MISTRAL PRO NC and MISTRAL PRO SC catheters have been designed, manufactured, tested and packaged with the utmost care. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Handling and storage of this device, as well as factors relating to the patient, diagnosis, treatment, surgical procedure and other areas beyond the control of HEXACATH, directly affect the device and the results obtained from its use. HEXACATH's obligations are limited to the replacement of the device. HEXACATH will not be liable in any way for any loss, damage or incidental or consequential costs arising from the use of the device. HEXACATH does not assume, nor does it authorise any third party to assume on its behalf, any additional liability or obligation in connection with the device.

HEXACATH shall not be liable for any reuse, reprocessing or re-sterilisation of the products and makes no warranty, express or implied, including but not limited to any warranty of merchantability or fitness for a particular purpose, with respect to these products.

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HEXACATH
4, passage St-Antoine
92500 Rueil-Malmaison - FRANCE
Tel. : +33 (0)1 41 39 01 92
Fax : +33 (0)1 41 39 04 13
www.hexacath.com

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² David E. Kandzari et al., "Clinical Outcomes Following Predilatation with a Novel 1.25-Mm Diameter Angioplasty Catheter", *Catheterization and Cardiovascular Interventions: Official Journal of the Society for Cardiac Angiography & Interventions* 77, n.º 4 (1 mars 2011): 510-14, <https://doi.org/10.1002/ccd.22734>

³ Erdem Özel et al., "What Is Better for Predilatation in Bioreabsorbable Vascular Scaffold Implantation: A Non-Compliant or a Compliant Balloon?", *Anatolian Journal of Cardiology* 16, n.º 4 (avril 2016): 244-49, <https://doi.org/10.5152/AnatJCardiol.2015.6184>

⁴ https://www.bostonscientific.com/content/dam/Manuals/us/current-rev-en/90990229-01A_Emerg_dfu_en-US_S.pdf

⁵ "Evaluation of Coronary Luminal Diameter Enlargement With Emerg™ 1.20 mm PTCA Dilatation Catheter - Study Results - ClinicalTrials.gov", <https://clinicaltrials.gov/ct2/show/results/NCT01635881>

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¹ "Sprinter Legend Balloon Catheter - Study Results ClinicalTrials.gov", <https://clinicaltrials.gov/ct2/show/results/NCT00961311>.

² David E. Kandzan et al., "Clinical Outcomes Following Predilatation with a Novel 1.25-Mm Diameter Angioplasty Catheter", *Catheterization and Cardiovascular Interventions: Official Journal of the Society for Cardiac Angiography & Interventions* 77, nr 4 (1 mars 2011): 510-14, <https://doi.org/10.11552/IntJCardiol.2015.6184>

³ Erdem Özel et al., "What Is Better for Predilatation in Bioreabsorbable Vascular Scaffold Implantation: A Non-Compliant or a Compliant Balloon?", *Anatolian Journal of Cardiology* 18, nr 4 (avril 2016): 244-49, <https://doi.org/10.5152/AnatolJCardiol.2015.6184>

⁴ https://www.bostonscientific.com/content/dam/Manuals/us/current-rev-en/90990229-01A_Emerge_dtu_en-US_S.pdf

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¹ "Sprinter Legend Balloon Catheter - Study ResultsClinicalTrials.gov", <https://clinicaltrials.gov/ct2/show/results/NCT00961311>.

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³ Erdem Özel et al., "What Is Better for Predilatation in Bioreabsorbable Vascular Scaffold Implantation: A Non-Compliant or a Compliant Balloon?", *Anatolian Journal of Cardiology* 16, nr 4 (avril 2016): 244-49, <https://doi.org/10.5152/AnatolJCardiol.2015.6184>

⁴ https://www.bostonscientific.com/content/dam/Manuals/us/current-rev-en/90990229-01A_Emerge_dtu_en-US_S.pdf

⁵ "Evaluation of Coronary Luminal Diameter Enlargement With Emerge™ 1.20 mm PTCA Dilatation Catheter - Study Results - ClinicalTrials.gov", <https://clinicaltrials.gov/ct2/show/results/NCT01635881>

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¹ "Sprinter Legend Balloon Catheter - Study ResultsClinicalTrials.Gov", <https://clinicaltrials.gov/ct2/show/results/NCT00961311>.

² David E. Kandzari et al., "Clinical Outcomes Following Predilatation with a Novel 1.25-Mm Diameter Angioplasty Catheter", *Catheterization and Cardiovascular Interventions: Official Journal of the Society for Cardiac Angiography & Interventions* 77, nr 4 (1 mars 2011): 510-14, <https://doi.org/10.1002/ccd.22734>

³ Erdem Özel et al., "What Is Better for Predilatation in Bioresorbable Vascular Scaffold Implantation: A Non-Compliant or a Compliant Balloon?", *Anatolian Journal of Cardiology* 16, nr 4 (avril 2016): 244-49, <https://doi.org/10.5152/AnatolJCardiol.2015.6184>

⁴ https://www.bostonscientific.com/content/dam/Manuals/us/current/rev-en/90990229-01A_Emerge_dfu_en-US_S.pdf

⁵ "Evaluation of Coronary Luminal Diameter Enlargement With Emerge™ 1.20 mm PTCA Dilatation Catheter - Study Results - ClinicalTrials.Gov", <https://clinicaltrials.gov/ct2/show/results/NCT01635881>

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² David E. Kandzari et al., "Clinical Outcomes Following Predilation with a Novel 1.25-Mm Diameter Angioplasty Catheter", *Catheterization and Cardiovascular Interventions: Official Journal of the Society for Cardiac Angiography & Interventions* 77, n° 4 (1 mars 2011): 510-14, <https://doi.org/10.1002/ccd.22734>

³ Erdem Özel et al., "What Is Better for Predilation in Bioresorbable Vascular Scaffold Implantation: A Non-Compliant or a Compliant Balloon?", *Anatolian Journal of Cardiology* 16, n° 4 (avril 2016): 244-49, <https://doi.org/10.5152/AnatolJCardiol.2015.6184>

⁴ https://www.bostonscientific.com/content/dam/Manuals/us/current-rev-en/90990229-01A_Emerge_dtu_en-US_S.pdf

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